



MODEL 190

Cryo-Can

Eliminates sample contamination during scanning electron microscope (SEM) operation.

Model 190 Cryo-Can Specifications

Capacity	270 ml liquid nitrogen
Hold time	Up to 6 hours
SEM operating conditions	Up to 30 keV
Flange interface	Fits all SEMs; subject to the availability of ports
Liquid Nitrogen vessel	Can be baked to 150 °C
SEM compatibility	Contact a Fischione Instruments representative to discuss SEM compatibility
Warranty	One year



E.A. Fischione Instruments, Inc.
9003 Corporate Circle
Export, PA 15632 USA
Tel: +1 724.325.5444
Fax: +1 724.325.5443
info@fischione.com
www.fischione.com

©2013 E.A. Fischione Instruments, Inc. All rights reserved.
Document Number SP190 Revision 00 07/2013