



MODEL 2040

Dual-Axis Tomography Holder

A holder that features an optimal tilt angle range in narrow gap (~ 5 mm) pole-piece geometries, while maintaining microscope resolution.

Model 2040 Dual-Axis Tomography Holder Specifications

Holder type	Room-temperature, single tilt, in-plane rotation For most transmission electron microscopes (TEMs)
Specimen size	3 mm diameter Wide range of thicknesses up to 100 μm
Resolution	0.34 nm (in all directions)
Drift	< 1.5 nm/minute
Maximum tilt range	Up to $\pm 70^\circ$
Field of view	Up to 950 μm at 70°
Warranty	One year

All specifications depend on the microscope model, pole piece type, and aperture position.
For ultimate resolution and drift performance, the TEM must meet the manufacturer's specifications.



E.A. Fischione Instruments, Inc.
9003 Corporate Circle
Export, PA 15632 USA
Tel: +1 724.325.5444
Fax: +1 724.325.5443
info@fischione.com
www.fischione.com

©2015 E.A. Fischione Instruments, Inc. All rights reserved.
Document Number SP2040 Revision 02 08/2015